


Fragment screening using Microscale Thermophoresis and X-ray Crystallography against the lysine methyltransferase G9a

domainex
the faster route to drug discovery

Introduction

- G9a is a lysine methyltransferase (KMT) involved in epigenetic gene regulation by covalent modification of histones
- G9a catalyses the transfer of methyl groups from S-adenosyl methionine (SAM) to lysine residues on histone proteins (Fig A)
- Literature supports the role of G9a in mechanisms of carcinogenesis, making it an attractive oncology target^[1-4]
- Domainex has solved the key technical drug discovery challenges associated with KMTs, including generating a number of proprietary crystal structures, assays and a novel screening library of small molecule inhibitors
- In this poster, we report a fragment-based hit-finding approach using Microscale Thermophoresis (MST) for the ternary G9a-SAM-fragment system and hit confirmation using protein X-ray crystallography



Microscale Thermophoresis (MST)

- MST detects changes in the size, charge and hydration shell of molecules and enables the measurement of biomolecule interactions or enzyme kinetics
- Fluorescently-labelled molecules are exposed to precise temperature gradients generated by an IR-laser within microliter-volume glass capillaries
- The thermophoretic properties of a molecule change upon ligand binding
- These changes can be used to derive dissociation constants (K_d) within minutes


Fragment Library

- Our fragment library contains ~1,000 fragments and has been carefully designed to maximise the chances of finding suitable chemical starting points
- The library gives good coverage of 'ideal fragment space' by optimising a number of parameters, is 'Rule of 3' compliant, has good chemical diversity and is SP3 rich
- Additional filters were applied to remove:
 - compounds containing atoms other than H, C, N, O, S, F, Cl and reactive functional groups
- All compounds in the library show > 1 mM aqueous solubility in 1% DMSO


G9a Fragment Screen

- We screened our fragment library at 1 mM against a G9a-SAM complex using MST. Fragments were declared as hits if a significant shift was observed (Fig B). The thermophoresis traces allow easy identification of false-positive fragments, which caused a shift in amplitude due to protein aggregation (Figs C and D, respectively).
- We obtained a 5.3% hit rate. Screening the same fragments using Differential Scanning Fluorimetry (DSF) or the activity-based AlphaScreen both showed only a 0.3% hit rate.
- Hits were taken into secondary screening to determine their binding affinities (K_d) to the G9a-SAM complex, again using MST (for example, see Fig. E)
- K_d values for 7 hits were also determined without the co-factor, SAM. The K_d s for two fragments were essentially unchanged whereas five fragments showed significant reduction in binding, highlighting the importance of being able to study a ternary system (Table F)
- Orthogonal confirmation of hit binding to G9a was demonstrated by Saturation Transfer Difference (STD) NMR spectroscopy (Fig. G)
- Three G9a-fragment structures were solved in-house in the presence of co-factor SAM with a resolution of 1.5- 2.0 Å (Fig H) which revealed different fragment binding modes
- This has led to several options for FBDD to provide alternative inhibitor chemotypes.


A) Schematic of G9a function




E) Dose-response measurements




B) MST single shot fragment screening, MTP#s given are hits



C) Thermophoresis traces




D) Aggregation



F) Screening summary


Frag ID	K_d + SAM	LE (+SAM)	K_d - SAM	Comment	STD-NMR	Crystal trials	X-ray Structure
	K_d [μ M]		K_d [μ M]	MOA	Positive binding	X-ray Structure	Resolution
MTP4E1	117	0.41	94	SAM independent	✓	X	
MTP3G1	718	0.36	518	SAM independent	X		
MTP3B6	17	0.65	109	SAM dependent	✓	✓	1.5 Å
MTP2C3	56	0.41	>1 mM	SAM dependent	X		
MTP3G10	195	0.56	>1 mM	SAM dependent	✓	✓	1.8 Å
MTP2D8	534	0.50	Non Binder	SAM dependent	✓	✓	2.0 Å
MTP2H9	564	0.44	Non Binder	SAM dependent	✓	X	

G) STD NMR confirmation



STD-NMR confirmation spectra for five fragments. Spectrum A (black) shows the NMR-STD confirmation for each fragment, as well as the SAM co-binder. Spectrum B (red) shows the reference spectrum of each fragment and the SAM co-binder (8.16 ppm) in PBS, pH 8.5 (10 % D₂O). Spectrum C (grey) is the false positive control. No signal indicates there was no aggregation of the fragment, and no direct excitation of the fragment or SAM with the on-resonance pulse.

H) In-house G9a-Fragment structures (Orange - Fragment, Yellow - SAM)



Summary

- We have used MST to successfully screen a sub-set of our fragment library against the KMT G9a with a hit rate of 5.3%, identifying fragments with high ligand efficiencies.
- MST allows easy false-positive identification by highlighting compounds that induce protein aggregation or cause fluorescent effects.
- We were able to identify low and high affinity binders using the same technique.
- Three fragment hits were successfully crystallised bound to G9a, which enabled a SBDD program for this target.

Services/Contact

Domainex welcomes interest from any potential collaborators, industrial or academic. If you would like to learn more about applying our drug-discovery platform to other targets, please contact: tom.mander@domainex.co.uk
www.domainex.co.uk

References

[1] Copeland et al., *Nature Reviews*, 2012, (8), 724-732; [2] Hamamoto et al., *Nature Cell Biology*, 2004, (6), 731-740; [3] Hamamoto et al., *Cancer Sci.*, 2006, (97), 113-118; [4] Liu et al. *J. Natl. Cancer. Inst.*, 2013, doi: 10.1093/jnci/djt30; [5] Pons et al. *Eur Heart J* 2009 Feb;30(3):266-77; [6] NanoTemper User Manual Monolith NT.115 2011 July; [7] Casciello et al., *Front Immunol.*, 2015, (6), 1-12