


PROTAC Platform: Enabling Rapid Design, Characterisation and Assessment

Alicia Galván Álvarez, Gabriel Negoita-Giras, Adam Wallington, Anna Hopkins, Aidan Johnson, Luke Williams, Megan O'Brian, Cameron Haddow, Sarah Whipple, Austra Jablonskyte, Natasha Pretot, Andrew Ratcliffe, Philip Fallon, Gary Pitt, Trevor Askwith and Catherine Stace, David Gibson, Mehul Jesani, Joe Mason, Sophie Williams

Domainex Ltd, Chesterford Research Park, Saffron Walden, CB10 1XL, UK; Domainex Ltd, Iconix 2, Unity Campus, London Road, Pampisford, CB22 3EG, UK

Introduction & Significance



What are PROTACs?

- Proteolysis targeting chimeras (**PROTACs**) are heterobifunctional molecules that induce proximity between an E3 ligase and protein of interest (POI) to facilitate ubiquitination and degradation of the target protein
- Ternary complex formation is required to induce proximity-driven ubiquitination
- Degradation of targeted protein is achieved by hijacking the ubiquitin-proteasome system which can offer numerous advantages versus traditional inhibitors

Significance of this study

- ✓ High-throughput parallel synthesis enables rapid cycle times and faster SAR
- ✓ An integrated platform promotes holistic understanding of chemical matter and aids smarter compound design in an area where property prediction is challenging
- ✓ Quick coverage of chemical space, establishment of SAR and hit ID, is better suited to earlier stage projects than traditional PROTAC platforms

Traditional PROTAC Platform

Phase I

Curation of POI ligand
Structural information
POI ligand SAR

Set of 10-20 designed PROTACs synthesised

- Different POI ligands
- Linker length & attachment group
- Choice of E3 ligases

ADME
Physchem
Profiling

Assays

- Hepatocyte turnover
- Biorelevant solubility
- Media cell stability
- Permeability
- Cell accumulation

Phase II

Optimise "best" PROTAC

Further SAR exploration around linker and linker attachments

- Potency
- Selectivity
- ADME

Nomination of PROTAC for *in vivo* studies

Biological
Profiling

Inhibition assays

- Degradation assays**
- Nano-Glo® HiBiT detection system
 - Automated western blotting

Evaluation of binary and ternary complex formation


- MST
- GCI
- FP
- Alpha screen/HTRF
- Nano DSF

Direct-to-Biology (D2B) PROTAC Platform


Plate-based synthesis of PROTACs




Automated UPLC analysis



Evaluation in degradation assays and hit ID



Hit Validation



Library of linker-E3 ligase ligands prepared in tube racks

Day 1: Set up reactions in 96 or 384 well plates using automation

Day 2: Dilute to assay concentration, generate QC plate

Day 7-8: Screening results available

Day 4-7: Cell based screening assay (Hit ID)

Day 3: Automated UPLC analysis, plates shipped for hit identification


Typical workflow and timelines for D2B approach

- Up to 384 compounds synthesised, analysed and screened in single plate format. Throughput of 1000 compounds/month
- Minimal reagent and solvent use (0.25 µmol POI ligand, 10 µL solvent/ reaction)
- Crude mixtures used in hit identification using cell-based Nano-Glo HiBiT assays
- Hit validation performed using automated western blotting system
- Validated hits to be resynthesized and purified for evaluation of binary/ternary complex formation and ADME

Biological Profiling


Degradation assays:

- Nano-Glo HiBiT assays can be used to quantify cellular protein and assess PROTAC potency (A)
- Native POI-specific protein degradation can be re-confirmed in an orthologous screen, using an automated, capillary-based immuno-assay system (Simple Western™) (B). Endogenous-tagging of POI is not required



Biophysical evaluation of binary and ternary complex formation:

- Biophysical suite (MST, GCI, NanoDSF, FP) for development of binding assays
- Characterise ligand and degrader binding to degradation receptor (A), POIs (B), and ternary complex formation (C)



Advancing the Field of PROTACs

- ✓ High-throughput plate-based chemistry (1000 compounds per month) and direct-to-biology to accelerate cycle times and rapidly build SAR
- ✓ Computational chemistry to direct initial PROTAC design phase (e.g. growth vectors for linkers, novel E3 ligase binding motifs and attachment points)
- ✓ Bespoke linker-E3 ligase ligand toolbox with high diversity (linker type, length and rigidity, E3 ligase type and attachment points, etc)
- ✓ Comprehensive biochemical, biophysical, and cellular profiling, integrated with chemical design & synthesis
- ✓ Rapid ADME and physchem profiling to accelerate property optimization

Domainex welcomes interest from any potential collaborators, industrial or academic. If you would like to learn more about our drug-discovery platform, please contact: enquiries@domainex.co.uk